




Canadian Angus Hair Sample Instructions

1



Pull hair sample from the switch of the tail. Include 50–60 hairs ensuring that the hair root bulbs are attached.

2



Place hair roots on the hair card (*there is a circle in the middle of the hair card that says Place Hair Roots Here*).

3



Affix hair with tape at the bottom of the card (*there is a rectangle that says Place Tape Here*).

4



Trim off excess hair extending from the hair card.

5



Label the hair card clearly, carefully, and accurately. Store hair card in a safe dry place until you send it to the lab.

6



Submit test request to CAA office. Send CAA-generated paperwork and hair card to Neogen Canada.

Canadian Angus Association
292140 Wagon Wheel Blvd,
Rocky View County, AB T4A 0E2

Neogen Canada
7323 Roper Road NW
Edmonton, AB T6E 0W4